

## Getting Started

### Before using LabQuest 2

- 1 Remove the pull tab, which protects the battery, from the back of the unit.
- 2 Connect the power adapter.
- 3 Charge for at least 12 hours.
- 4 Read the safety information found in this Quick-Start Guide.

Complete user guide is available at  
[www.vernier.com/labquest/guide](http://www.vernier.com/labquest/guide)

### Quick Data Collection

#### 1 Turn on LabQuest 2.

#### 2 Connect a sensor.

The sensor auto-IDs and displays a live readout. The default collection rate for the sensor is set up.

#### 3 Tap Collect ➤.

You are now collecting data. LabQuest App switches to the Graph screen when data collection begins.

## Safety Information

	Read all safety information and operating instructions included in this Quick-Start Guide prior to using LabQuest 2.
	LabQuest 2 is designed to be splash resistant. However, avoid water immersion and standing liquid on the display. If water gets in the device, immediately shut down the device (tap System on the Home screen, then tap Shut Down, or hold down the Power button until the shutdown sequence begins). Remove the battery, connected cables, SD card, and any other accessories. Allow to dry thoroughly before restarting. Do not attempt to dry using an external heat source.
	Safe operating temperatures are from 0°C to 45°C. Storage temperatures are -30°C to 60°C. Exposures to low or high extreme temperatures will temporarily reduce battery life. Avoid rapid temperature changes as condensation may form inside the device. Do not leave in a car, as temperatures can exceed the maximum storage temperature.
	LabQuest 2 contains a lithium-ion battery. Use only the supplied battery for this device. Do not puncture or expose to excessive heat or flame.
	Do not store LabQuest 2 in a chemical closet or in areas of concentrated chemical gases.

## Care of LabQuest

**Screen Calibration** – The screen may need to be recalibrated periodically. To access the calibration tool, tap System on the Home screen, then tap Calibrate Screen. If you are unable to access the Calibration tool, press and hold the Home button to display the Calibration tool. Tap the calibration target with the stylus in the five locations indicated to calibrate the screen.

**System Reset** – To reset the system, tap System on the Home screen, then tap Reboot. If you are unable to access the Reboot tool, press and hold the Power button. When the LabQuest shutdown message is displayed, release the Power button. Wait until the system shuts down, then turn the unit back on.

**Battery Maintenance** – Use only the supplied AC adapter or optional Charging Station to charge the LabQuest battery. A full charge can take twelve hours. The battery cannot be overcharged, and it can be safely recharged after a partial discharge.

**Screen Maintenance** – The LabQuest screen is splash resistant. Wipe the screen clean with a cotton cloth that is slightly dampened with water or ethanol. Do not use any other solvents. Do not submerge LabQuest in liquids.

## FEDERAL COMMUNICATIONS COMMISSION INTERFERENCE STATEMENT

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

This device complies with Industry Canada license requirements.  
(1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

**Industry Canada - Class B** This digital apparatus complies with Canadian ICES-003.  
as set out in the interference-causing equipment notice issued by Industry Canada.  
Operation is subject to the following two conditions:  
any interference, including interference that may cause undesired operation of the device.

- Reorient or relocate the receiving antenna.
  - Increase the separation between the equipment and receiver.
  - Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
  - Consult the dealer or an experienced radio/TV technician for help.

#### ECC Caution

**FCC caution:**  
This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:  
(1) this device may not cause harmful interference and  
(2) this device must accept any interference received, including interference that may cause undesired operation

RF Exposure Warning

**RF Exposure Warning**  
The equipment complies with RF exposure limits set forth for an uncontrolled environment. The antenna(s) used for this transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

You are cautioned that changes or modifications not expressly approved by the party responsible for compliance could void your authority to operate the equipment.

not RSS standard(s). Operation is subject to the following two conditions:  
device must accept any interference, including interference that may cause

not exceed the Class B limits for radio noise emissions from digital apparatus as per Part 15 of the Canadian Electrical Code and entitled "Digital Apparatus," ICES-003 of Industry Canada.

(1) this device may not cause interference, and (2) this device must accept undesired operation of the device.

The antenna type and its gain should be so chosen that the equivalent exposure limits set forth for an uncontrolled environment. The antenna(s) operating in conjunction with any other antenna or transmitter

*Canada applicables aux appareils radio exempts de licence. L'exploitation est  
permise sous les deux conditions suivantes : (1) l'appareil ne doit pas produire de brouillage, et (2) l'appareil doit accepter tout  
brouillage susceptible d'en compromettre le fonctionnement.*

**électriques applicables aux appareils numériques de Classe B prescrites**  
**appareils Numériques," NMB-003 édictée par industrie Canada.**

*cet appareil ne peut causer d'interférences, et (2) cet appareil doit accepter  
toute interférence, y compris celles qui peuvent provoquer un fonctionnement du dispositif.*

es autres utilisateurs, le type d'antenne et son gain doivent être choisis de manière à ce que l'écart entre la puissance émise et la puissance reçue (E.i.r.p) n'est pas plus grand que celui permis pour une communication établie. Ces limites sont fixées par les normes de sécurité et doivent être respectées pour éviter les interférences avec d'autres appareils. Les appareils doivent également être placés de manière à ce qu'ils ne soient pas exposés à des sources de champs électriques ou magnétiques qui pourraient perturber leur fonctionnement.



 LABQUEST<sup>®</sup> 2 Quick

What

- LabC
  - Rech
  - Powe
  - USB
  - Stylu
  - Stylu

Logger List  
(optional)  
at [www.yourlogger.com](http://www.yourlogger.com)



[Warranty Info](#)

**LabQuest 2** – This product is warranted to be free from defects in materials and workmanship for a period of five years from the date of shipment.

**LabQuest 2 Battery** – This product is warranted to be free from defects in materials and workmanship for a period of one year from the date of shipment.

Warranty covers use by educational institution only. It does not cover damage to the product caused by abuse or improper use.